

## Biotinylated Anti-bovine IL-17A Polyclonal Antibody

KPAB0036

### Description

---

This Biotinylated Anti-bovine IL-17A Polyclonal Antibody is has been validated for ELISA, ELISpot applications.

### Product Information

---

<b>SKU:</b>	KPAB0036
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Biotin
<b>Applications:</b>	ELISA, ELISpot

### Antibody Data

---

<b>Host:</b>	Rabbit
<b>Reactivity:</b>	Bovine
<b>Immunogen:</b>	Recombinant Bovine IL-17A
<b>Molecular Weight:</b>	17.4 kDa (calculated) for bovine IL-17A
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN <sub>3</sub>
<b>Homology:</b>	Bos indicus (zebu), Bos javanicus (banteng), Bos taurus (cattle)
<b>Cross-Reactivity:</b>	Canine IL-17A: Weak; Caprine IL-17A: Moderate; Cynomolgus IL-17A: Weak; Dolphin IL-17A: Strong; Feline IL-17A: Moderate; Ferret IL-17A: Weak; Guinea Pig IL-17A: Weak; Human IL-17A: Weak; Mouse IL-17A: Weak; Ovine IL-17A: Strong; Rabbit IL-17A: Weak; Rat IL-17A: Weak; Swine IL-17A: Strong

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.