

Biotinylated Anti-bovine VEGF-A Polyclonal Antibody

KPAB0038

Description

This Biotinylated Anti-bovine VEGF-A Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0038
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Bovine
Immunogen:	Recombinant Bovine VEGF-A
Molecular Weight:	19.2 kDa (calculated) for bovine VEGF-A
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN3
Cross-Reactivity:	Canine VEGF-A: Strong; Equine VEGF-A: Strong; Feline VEGF-A: Strong; Mouse VEGF-A: Strong; Ovine VEGF-A: Strong; Rabbit VEGF-A: Weak; Rabbit VEGF-A165: Weak; Rat VEGF-A: Strong; Swine VEGF-A: Strong
Homology:	Bison bison bison (American buffalo), Bos taurus (cattle), Bubalus bubalis (water buffalo), Bubalus kerabau (carabao)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.