

Anti-feline IFN gamma Polyclonal Antibody

KPAB0040

Description

This Anti-feline IFN gamma Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0040
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Feline
Immunogen:	Recombinant Feline IFN gamma
Molecular Weight:	17.0 kDa (calculated) for feline IFN- η
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Canine IFN-H: Weak
Homology:	Acinonyx jubatus (cheetah), Felis catus (domestic cat), Prionailurus bengalensis (leopard cat), Prionailurus viverrinus (fishing cat), Puma concolor (puma)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.