

Biotinylated Anti-feline IFN gamma Polyclonal Antibody

KPAB0042

Description

This Biotinylated Anti-feline IFN gamma Polyclonal Antibody has been validated for WB, ELISA applications.

Product Information

SKU: KPAB0042

Conjugated: Biotin

Applications: WB, ELISA

Antibody Data

Host: Rabbit

Reactivity: Feline

Immunogen: Recombinant Feline IFN gamma

Purification: Antigen-affinity Purification

Homology: Acinonyx jubatus (cheetah), Felis catus (domestic cat), Prionailurus bengalensis (leopard cat), Prionailurus viverrinus (fishing cat), Puma concolor (puma)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.