

Anti-equine IL-17A Polyclonal Antibody

KPAB0045

Description

This Anti-equine IL-17A Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0045
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Equine
Immunogen:	Recombinant Horse IL-17A
Molecular Weight:	14.9 kDa (calculated) for equine IL-17A
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN3
Cross-Reactivity:	Bovine IL-17A: Strong; Swine IL-17A: Strong
Homology:	Equus caballus (horse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.