

Anti-swine IL-22 Polyclonal Antibody

KPAB0049

Description

This Anti-swine IL-22 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0049
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Swine
Immunogen:	Recombinant Porcine IL-22
Molecular Weight:	18.0 kDa (calculated) for swine IL-22
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Homology:	Phacochoerus africanus (Common warthog), Sus scrofa (pig)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.