

Biotinylated Anti-rabbit IL-4 Polyclonal Antibody

KPAB0054

Description

This Biotinylated Anti-rabbit IL-4 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0054
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Goat
Reactivity:	Rabbit
Immunogen:	Recombinant Rabbit IL-4
Molecular Weight:	14.2 kDa (calculated) for rabbit IL-4
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Homology:	Oryctolagus cuniculus (rabbit)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.