

Anti-swine IL-6 Polyclonal Antibody

KPAB0057

Description

This Anti-swine IL-6 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

| | |
|-----------------------|--|
| SKU: | KPAB0057 |
| Concentration: | 1.0 mg/mL |
| Conjugated: | Unconjugated |
| Applications: | ELISA, ELISpot, FC, Neutralization, WB |

Antibody Data

| | |
|--------------------------|---|
| Host: | Rabbit |
| Reactivity: | Swine |
| Immunogen: | Recombinant Porcine IL-6 |
| Molecular Weight: | 20.9 kDa (calculated) for swine IL-6 |
| Purification: | Antigen-affinity Purification |
| Formulation: | 1X PBS buffer containing 0.09% NaN ₃ |
| Cross-Reactivity: | Feline IL-6: Weak |
| Homology: | Sus scrofa (pig) |

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.