

Anti-bovine IL-1F5 Polyclonal Antibody

KPAB0063

Description

This Anti-bovine IL-1F5 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0063
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Bovine
Immunogen:	Recombinant Bovine IL-1F5
Molecular Weight:	17.1 kDa (calculated) for bovine IL-1F5
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Homology:	Bos indicus (zebu), Bos taurus (cattle)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.