

Equine IL-17A Monoclonal Antibody (clone K0078E-1)

KPAB0067

Description

This Equine IL-17A Monoclonal Antibody (clone K0078E-1) has been validated for ELISA, FC applications.

Product Information

SKU: KPAB0067
Conjugated: Unconjugated
Applications: ELISA, FC

Antibody Data

Host: Mouse
Reactivity: Equine
Immunogen: Recombinant Horse IL-17A
Purification: Dialysis
Homology: Equus caballus (horse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.