

## Equine IL-17A Monoclonal Antibody (clone K0078E-2)

**KPAB0069**

### Description

---

This Equine IL-17A Monoclonal Antibody (clone K0078E-2) has been validated for ELISA, FC applications.

### Product Information

---

**SKU:** KPAB0069  
**Conjugated:** Unconjugated  
**Applications:** ELISA, FC

### Antibody Data

---

**Host:** Mouse  
**Reactivity:** Equine  
**Immunogen:** Recombinant Horse IL-17A  
**Purification:** Dialysis  
**Homology:** Equus caballus (horse)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.