

Biotinylated Equine IL-17A Monoclonal Antibody (clone K0078E-1)

KPAB0070

Description

This Biotinylated Equine IL-17A Monoclonal Antibody (clone K0078E-1) has been validated for ELISA applications.

Product Information

SKU: KPAB0070

Conjugated: Biotin

Applications: ELISA

Antibody Data

Host: Mouse

Reactivity: Equine

Immunogen: Recombinant Horse IL-17A

Purification: Dialysis

Homology: Equus caballus (horse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.