

## Biotinylated Equine IL-17A Monoclonal Antibody (clone K0078E-1)

**KPAB0071**

### Description

---

This Biotinylated Equine IL-17A Monoclonal Antibody (clone K0078E-1) has been validated for ELISA applications.

### Product Information

---

**SKU:** KPAB0071

**Conjugated:** Biotin

**Applications:** ELISA

### Antibody Data

---

**Host:** Mouse

**Reactivity:** Equine

**Immunogen:** Recombinant Horse IL-17A

**Purification:** Dialysis

**Homology:** Equus caballus (horse)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.