

## Anti-rabbit IFN gamma Polyclonal Antibody

**KPAB0076**

### Description

---

This Anti-rabbit IFN gamma Polyclonal Antibody has been validated for ELISA applications.

### Product Information

---

**SKU:** KPAB0076  
**Conjugated:** Unconjugated  
**Applications:** ELISA

### Antibody Data

---

**Host:** Goat  
**Reactivity:** Rabbit  
**Immunogen:** Recombinant Rabbit IFN gamma  
**Purification:** Antigen-affinity Purification

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.