

Anti-equine VEGF-A Polyclonal Antibody

KPAB0080

Description

This Anti-equine VEGF-A Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0080
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Goat
Reactivity:	Equine
Immunogen:	Recombinant Horse VEGF-A
Molecular Weight:	19.2 kDa (calculated) for equine VEGF-A
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN3
Cross-Reactivity:	Bovine VEGF-A: Weak; Canine VEGF-A: Weak; Mouse VEGF-A: Weak; Rat VEGF-A: Weak; Swine VEGF-A: Moderate
Homology:	Equus asinus (ass), Equus caballus (horse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.