

## Anti-canine TNF alpha Polyclonal Antibody

KPAB0081

### Description

---

This Anti-canine TNF alpha Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

### Product Information

---

<b>SKU:</b>	KPAB0081
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Unconjugated
<b>Applications:</b>	ELISA, ELISpot, FC, Neutralization, WB

### Antibody Data

---

<b>Host:</b>	Goat
<b>Reactivity:</b>	Canine
<b>Immunogen:</b>	Recombinant Canine TNF alpha
<b>Molecular Weight:</b>	17.3 kDa (calculated) for canine TNF- $\alpha$
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN <sub>3</sub>
<b>Cross-Reactivity:</b>	Cynomolgus TNF-A: Strong; Equine TNF-A: Weak; Feline TNF-A: Moderate; Human TNF-A: Strong; Rabbit TNF-A: Moderate
<b>Homology:</b>	Canis lupus dingo (dingo), Canis lupus familiaris (dog), Nyctereutes procyonoides (raccoon dog)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.