

Biotinylated Anti-swine CCL3L1 Polyclonal Antibody

KPAB0091

Description

This Biotinylated Anti-swine CCL3L1 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0091
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Swine
Immunogen:	Recombinant Porcine CCL3L1
Molecular Weight:	7.8 kDa (calculated) for swine CCL3L1
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Equine CCL3: Moderate
Homology:	Sus scrofa (pig)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.