

Anti-equine CCL3 (MIP-1 alpha) Polyclonal Antibody

KPAB0092

Description

This Anti-equine CCL3 (MIP-1 alpha) Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0092
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Equine
Immunogen:	Recombinant Horse CCL3
Molecular Weight:	7.9 kDa (calculated) for equine CCL3
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Homology:	Equus caballus (horse), Equus przewalskii (Przewalski's horse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.