

## Equine IL-4 Monoclonal Antibody (clone K0003E-1)

**KPAB0094**

### Description

---

This Equine IL-4 Monoclonal Antibody (clone K0003E-1) has been validated for ELISA, FC applications.

### Product Information

---

**SKU:** KPAB0094

**Conjugated:** Unconjugated

**Applications:** ELISA, FC

### Antibody Data

---

**Host:** Mouse

**Reactivity:** Equine

**Immunogen:** Recombinant Horse IL-4

**Purification:** Dialysis

**Homology:** Equus asinus (ass), Equus caballus (horse), Equus przewalskii (Przewalski's horse), Equus zebra hartmannae (Hartmann's mountain zebra)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.