

Equine IL-4 Monoclonal Antibody (clone K0003E-1)

KPAB0095

Description

This Equine IL-4 Monoclonal Antibody (clone K0003E-1) has been validated for ELISA, FC applications.

Product Information

SKU: KPAB0095
Conjugated: Unconjugated
Applications: ELISA, FC

Antibody Data

Host: Mouse
Reactivity: Equine
Immunogen: Recombinant Horse IL-4
Purification: Dialysis
Homology: Equus asinus (ass), Equus caballus (horse), Equus przewalskii (Przewalski's horse), Equus zebra hartmannae (Hartmann's mountain zebra)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.