

## Anti-feline TNF alpha Polyclonal Antibody

### KPAB0103

#### Description

---

This Anti-feline TNF alpha Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

#### Product Information

---

<b>SKU:</b>	KPAB0103
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Unconjugated
<b>Applications:</b>	ELISA, ELISpot, FC, Neutralization, WB

#### Antibody Data

---

<b>Host:</b>	Rabbit
<b>Reactivity:</b>	Feline
<b>Immunogen:</b>	Recombinant Feline TNF alpha
<b>Molecular Weight:</b>	17.2 kDa (calculated) for feline TNF- $\alpha$
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN <sub>3</sub>
<b>Cross-Reactivity:</b>	Canine TNF-A: Strong; Caprine TNF-A: Weak; Cynomolgus Monkey TNF-A: Strong; Equine TNF-A: Weak; Human TNF-A: Moderate; Swine TNF-A: Weak
<b>Homology:</b>	Felis catus (domestic cat), Prionailurus bengalensis (leopard cat), Prionailurus viverrinus (fishing cat)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.