

## Biotinylated Anti-equine CXCL10 (IP-10) Polyclonal Antibody

**KPAB0107**

### Description

---

This Biotinylated Anti-equine CXCL10 (IP-10) Polyclonal Antibody has been validated for ELISA, ELISpot applications.

### Product Information

---

<b>SKU:</b>	KPAB0107
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Biotin
<b>Applications:</b>	ELISA, ELISpot

### Antibody Data

---

<b>Host:</b>	Rabbit
<b>Reactivity:</b>	Equine
<b>Immunogen:</b>	Recombinant Horse CXCL10
<b>Molecular Weight:</b>	9.4 kDa (calculated) for equine CXCL10
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN3
<b>Cross-Reactivity:</b>	Bovine CXCL10: Weak; Dolphin CXCL10: Weak; Ovine CXCL10: Weak; Swine CXCL10: Weak
<b>Homology:</b>	Equus caballus (horse), Equus przewalskii (Przewalski's horse)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.