

Anti-dolphin TNF alpha Polyclonal Antibody

KPAB0114

Description

This Anti-dolphin TNF alpha Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0114
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Goat
Reactivity:	Dolphin
Immunogen:	Recombinant Dolphin TNF alpha
Molecular Weight:	17.4 kDa (calculated) for dolphin TNF- α
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Bovine TNF-A: Weak; Rabbit TNF-A: Weak; Swine TNF-A: Weak
Homology:	Orcinus orca (killer whale), Steno bredanensis (rough-toothed dolphin), Tursiops truncatus (common bottlenose dolphin)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.