

Anti-swine IFN alpha 1 Polyclonal Antibody

KPAB0116

Description

This Anti-swine IFN alpha 1 Polyclonal Antibody has been validated for WB, ELISA applications.

Product Information

SKU: KPAB0116
Conjugated: Unconjugated
Applications: WB, ELISA

Antibody Data

Host: Goat
Reactivity: Swine
Immunogen: Recombinant Porcine IFN alpha 1
Purification: Antigen-affinity Purification

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.