

Biotinylated Anti-dolphin TNF alpha Polyclonal Antibody

KPAB0120

Description

This Biotinylated Anti-dolphin TNF alpha Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0120
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Goat
Reactivity:	Dolphin
Immunogen:	Recombinant Dolphin TNF alpha
Purification:	Antigen-affinity Purification
Homology:	Orcinus orca (killer whale), Steno bredanensis (rough-toothed dolphin), Tursiops truncatus (common bottlenose dolphin)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.