

Biotinylated Anti-swine IFN alpha 1 Polyclonal Antibody

KPAB0122

Description

This Biotinylated Anti-swine IFN alpha 1 Polyclonal Antibody has been validated for ELISA applications.

Product Information

SKU: KPAB0122

Conjugated: Biotin

Applications: ELISA

Antibody Data

Host: Goat

Reactivity: Swine

Immunogen: Recombinant Porcine IFN alpha 1

Purification: Antigen-affinity Purification

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.