

Biotinylated Anti-bovine IL-15 Polyclonal Antibody

KPAB0124

Description

This Biotinylated Anti-bovine IL-15 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0124
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Goat
Reactivity: Bovine
Immunogen: Recombinant Bovine IL-15
Purification: Antigen-affinity Purification
Homology: Bos indicus (zebu), Bos mutus (wild yak), Bos taurus (cattle), Bubalus bubalis (water buffalo), Bubalus kerabau (carabao), Moschus berezovskii (Chinese forest musk deer)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.