

Equine IFN gamma Monoclonal Antibody (clone K0012E-1)

KPAB0126

Description

This Equine IFN gamma Monoclonal Antibody (clone K0012E-1) has been validated for ELISA, FC applications.

Product Information

SKU:	KPAB0126
Conjugated:	Unconjugated
Applications:	ELISA, FC

Antibody Data

Host:	Mouse
Reactivity:	Equine
Immunogen:	Recombinant Horse IFN gamma
Purification:	Dialysis
Homology:	Equus caballus (horse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.