

## Biotinylated Equine IFN gamma Monoclonal Antibody (clone K0012E-1)

### KPAB0128

#### Description

---

This Biotinylated Equine IFN gamma Monoclonal Antibody (clone K0012E-1) has been validated for ELISA applications.

#### Product Information

---

**SKU:** KPAB0128

**Conjugated:** Biotin

**Applications:** ELISA

#### Antibody Data

---

**Host:** Mouse

**Reactivity:** Equine

**Immunogen:** Recombinant Horse IFN gamma

**Purification:** Dialysis

**Homology:** Equus caballus (horse)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.