

## Anti-bovine CCL4 (MIP-1 beta) Polyclonal Antibody

**KPAB0130**

### Description

---

This Anti-bovine CCL4 (MIP-1 beta) Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

### Product Information

---

<b>SKU:</b>	KPAB0130
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Unconjugated
<b>Applications:</b>	ELISA, ELISpot, FC, Neutralization, WB

### Antibody Data

---

<b>Host:</b>	Goat
<b>Reactivity:</b>	Bovine
<b>Immunogen:</b>	Recombinant Bovine CCL4
<b>Molecular Weight:</b>	7.8 kDa (calculated) for bovine CCL4
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN <sub>3</sub>
<b>Cross-Reactivity:</b>	Canine CCL4: Strong; Dolphin CCL4: Weak; Equine CCL4: Moderate; Human CCL4: Moderate; Rabbit CCL4: Strong; Rat CCL4: Strong; Swine CCL4: Strong
<b>Homology:</b>	Bison bison bison (American buffalo), Bos indicus (zebu), Bos mutus (wild yak), Bos taurus (cattle)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.