

Anti-rat APRIL (TNFSF13) Polyclonal Antibody

KPAB0132

Description

This Anti-rat APRIL (TNFSF13) Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0132
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Goat
Reactivity:	Rat
Immunogen:	Recombinant Rat APRIL
Molecular Weight:	16.4 kDa (calculated) for rat APRIL
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN3
Cross-Reactivity:	Bovine APRIL: Weak; Mouse APRIL: Strong
Homology:	Rattus norvegicus (Norway rat)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.