

Anti-bovine IFN alpha A Polyclonal Antibody

KPAB0134

Description

This Anti-bovine IFN alpha A Polyclonal Antibody has been validated for WB, ELISA applications.

Product Information

SKU: KPAB0134
Conjugated: Unconjugated
Applications: WB, ELISA

Antibody Data

Host: Goat
Reactivity: Bovine
Immunogen: Recombinant Bovine IFN alpha A
Purification: Antigen-affinity Purification

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.