

## Biotinylated Anti-bovine CCL4 (MIP-1 beta) Polyclonal Antibody

**KPAB0138**

### Description

---

This Biotinylated Anti-bovine CCL4 (MIP-1 beta) Polyclonal Antibody has been validated for ELISA, ELISpot applications.

### Product Information

---

**SKU:** KPAB0138

**Conjugated:** Biotin

**Applications:** ELISA, ELISpot

### Antibody Data

---

**Host:** Goat

**Reactivity:** Bovine

**Immunogen:** Recombinant Bovine CCL4

**Purification:** Antigen-affinity Purification

**Homology:** Bison bison bison (American buffalo), Bos indicus (zebu), Bos mutus (wild yak), Bos taurus (cattle)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.