

Biotinylated Anti-bovine IFN alpha A Polyclonal Antibody

KPAB0142

Description

This Biotinylated Anti-bovine IFN alpha A Polyclonal Antibody has been validated for ELISA applications.

Product Information

SKU: KPAB0142

Conjugated: Biotin

Applications: ELISA

Antibody Data

Host: Goat

Reactivity: Bovine

Immunogen: Recombinant Bovine IFN alpha A

Purification: Antigen-affinity Purification

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.