

Anti-mouse IL-17A Polyclonal Antibody

KPAB0147

Description

This Anti-mouse IL-17A Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0147
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Goat
Reactivity:	Mouse
Immunogen:	Recombinant Mouse IL-17A
Molecular Weight:	15.0 kDa (calculated) for mouse IL-17A
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Equine IL-17A: Weak; Rabbit IL-17A: Weak; Rat IL-17A: Strong
Homology:	Mus musculus (house mouse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.