

Biotinylated Anti-mouse TNF alpha Polyclonal Antibody

KPAB0158

Description

This Biotinylated Anti-mouse TNF alpha Polyclonal Antibody has been validated for ELISA applications.

Product Information

SKU: KPAB0158

Conjugated: Biotin

Applications: ELISA

Antibody Data

Host: Rabbit

Reactivity: Mouse

Immunogen: Recombinant Mouse TNF alpha

Purification: Antigen-affinity Purification

Homology: Mus musculus (house mouse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.