

Anti-mouse VEGF-A Polyclonal Antibody

KPAB0159

Description

This Anti-mouse VEGF-A Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0159
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Mouse
Immunogen:	Recombinant Mouse VEGF-A
Molecular Weight:	17.3 kDa (calculated) for mouse TNF- α
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.02% NaN ₃
Homology:	Mus musculus (house mouse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.