

Botinylated Anti-mouse VEGF-A Polyclonal Antibody

KPAB0160

Description

This Botinylated Anti-mouse VEGF-A Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0160
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Mouse
Immunogen:	Recombinant Mouse VEGF-A
Molecular Weight:	19.3 kDa (calculated) for mouse VEGF-A
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.02% NaN3
Cross-Reactivity:	Bovine VEGF-A: Strong; Canine VEGF-A: Strong; Caprine VEGF-A: Strong; Chicken/Turkey VEGF-A: Weak; Equine VEGF-A: Strong; Feline VEGF-A: Strong; Human/Monkey VEGF-A: Strong; Ovine VEGF-A: Strong; Rat VEGF-A: Strong; Swine VEGF-A: Strong
Homology:	Mus musculus (house mouse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.