

Anti-canine IL-1 Receptor Antagonist Polyclonal Antibody

KPAB0162

Description

This Anti-canine IL-1 Receptor Antagonist Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0162
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Goat
Reactivity:	Canine
Immunogen:	Recombinant Canine IL-1 Receptor Antagonist
Molecular Weight:	16.7 kDa (calculated) for canine IL-1RA
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Dolphin IL-1Ra: Weak; Equine IL-1Ra: Weak; Mouse IL-1Ra: Moderate; Swine IL-1Ra: Weak
Homology:	Canis lupus dingo (dingo), Canis lupus familiaris (dog), Nyctereutes procyonoides (raccoon dog), Vulpes lagopus (Arctic fox), Vulpes vulpes (red fox)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.