

## Anti-ovine VEGF-A Polyclonal Antibody

KPAB0171

### Description

---

This Anti-ovine VEGF-A Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

### Product Information

---

<b>SKU:</b>	KPAB0171
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Unconjugated
<b>Applications:</b>	ELISA, ELISpot, FC, Neutralization, WB

### Antibody Data

---

<b>Host:</b>	Goat
<b>Reactivity:</b>	Ovine
<b>Immunogen:</b>	Recombinant Sheep VEGF-A
<b>Molecular Weight:</b>	14.1 kDa (calculated) for ovine VEGF-A
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN <sub>3</sub>
<b>Cross-Reactivity:</b>	Bovine VEGF-A: Moderate; Canine VEGF-A: Weak; Equine VEGF-A: Weak; Feline VEGF-A: Weak; Mouse VEGF-A: Weak; Rat VEGF-A: Weak; Swine VEGF-A: Moderate
<b>Homology:</b>	Oryx dammah (scimitar-horned oryx), Ovis aries (sheep), Pantholops hodgsonii (chiru)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.