

Biotinylated Anti-ovine VEGF-A Polyclonal Antibody

KPAB0172

Description

This Biotinylated Anti-ovine VEGF-A Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

| | |
|-----------------------|----------------|
| SKU: | KPAB0172 |
| Concentration: | 1.0 mg/mL |
| Conjugated: | Biotin |
| Applications: | ELISA, ELISpot |

Antibody Data

| | |
|--------------------------|--|
| Host: | Goat |
| Reactivity: | Ovine |
| Immunogen: | Recombinant Sheep VEGF-A |
| Molecular Weight: | 14.1 kDa (calculated) for ovine VEGF-A |
| Purification: | Antigen-affinity Purification |
| Formulation: | 1X PBS buffer containing 0.09% NaN ₃ |
| Cross-Reactivity: | Bovine VEGF-A: Moderate; Canine VEGF-A: Weak; Equine VEGF-A: Weak; Feline VEGF-A: Weak; Mouse VEGF-A: Weak; Rat VEGF-A: Weak; Swine VEGF-A: Moderate |
| Homology: | Oryx dammah (scimitar-horned oryx), Ovis aries (sheep), Pantholops hodgsonii (chiru) |

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.