

Biotinylated Anti-rabbit CCL2 (MCP-1) Polyclonal Antibody

KPAB0174

Description

This Biotinylated Anti-rabbit CCL2 (MCP-1) Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0174
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Goat
Reactivity:	Rabbit
Immunogen:	Recombinant Rabbit CCL2
Molecular Weight:	11.4 kDa (calculated) for rabbit CCL2
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Bovine CCL2: Weak; Canine CCL2: Weak; Equine CCL2: Weak; Mouse CCL2: Weak; Swine CCL2: Weak
Homology:	Oryctolagus cuniculus (rabbit)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.