

Anti-rabbit CXCL10 (IP-10) Polyclonal Antibody

KPAB0175

Description

This Anti-rabbit CXCL10 (IP-10) Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0175
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Goat
Reactivity:	Rabbit
Immunogen:	Recombinant Rabbit CXCL10
Molecular Weight:	8.6 kDa (calculated) for rabbit CXCL10
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Canine CXCL10: Weak; Dolphin CXCL10: Weak; Equine CXCL10: Weak; Feline CXCL10: Weak; Human CXCL10: Weak; Mouse CXCL10: Weak; Swine CXCL10: Weak
Homology:	Oryctolagus cuniculus (rabbit)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.