

Biotinylated Anti-mouse IL-22 Polyclonal Antibody

KPAB0180

Description

This Biotinylated Anti-mouse IL-22 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0180
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Chicken
Reactivity:	Mouse
Immunogen:	Recombinant Mouse IL-22
Molecular Weight:	16.7 kDa (calculated) for mouse IL-22
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Bovine IL-22: Moderate
Homology:	Mus musculus (house mouse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.