

Anti-bovine CCL11 (Eotaxin) Polyclonal Antibody

KPAB0181

Description

This Anti-bovine CCL11 (Eotaxin) Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0181
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Chicken
Reactivity:	Bovine
Immunogen:	Recombinant Bovine CCL11
Molecular Weight:	8.6 kDa (calculated) for bovine CCL11
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Equine CCL11: Weak; Guinea Pig CCL11: Weak
Homology:	Bison bison bison (American buffalo), Bos indicus (zebu), Bos javanicus (banteng), Bos mutus (wild yak), Bos taurus (cattle), Bubalus bubalis (water buffalo), Bubalus kerabau (carabao)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.