

Biotinylated Anti-dolphin IL-1 Receptor Antagonist Polyclonal Antibody

KPAB0184

Description

This Biotinylated Anti-dolphin IL-1 Receptor Antagonist Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

| | |
|-----------------------|----------------|
| SKU: | KPAB0184 |
| Concentration: | 1.0 mg/mL |
| Conjugated: | Biotin |
| Applications: | ELISA, ELISpot |

Antibody Data

| | |
|--------------------------|---|
| Host: | Rabbit |
| Reactivity: | Dolphin |
| Immunogen: | Recombinant Dolphin IL-1 Receptor Antagonist |
| Molecular Weight: | 17.0 kDa (calculated) for dolphin IL-1RA |
| Purification: | Antigen-affinity Purification |
| Formulation: | 1X PBS buffer containing 0.09% NaN ₃ |
| Cross-Reactivity: | Swine IL-1Ra: Weak |
| Homology: | Delphinapterus leucas (beluga whale), Globicephala melas (long-finned pilot whale), Monodon monoceros (narwhal), Phocoena sinus (vaquita), Tursiops truncatus (common bottlenose dolphin) |

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.