

## Biotinylated Anti-equine IFN beta 1 Polyclonal Antibody

### KPAB0192

#### Description

---

This Biotinylated Anti-equine IFN beta 1 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

#### Product Information

---

**SKU:** KPAB0192  
**Conjugated:** Biotin  
**Applications:** ELISA, ELISpot

#### Antibody Data

---

**Host:** Rabbit  
**Reactivity:** Equine  
**Immunogen:** Recombinant Horse IFN beta 1  
**Purification:** Antigen-affinity Purification  
**Homology:** Equus caballus (horse)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.