

Biotinylated Anti-swine Erythropoietin Polyclonal Antibody

KPAB0194

Description

This Biotinylated Anti-swine Erythropoietin Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0194
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Swine
Immunogen: Recombinant Porcine Erythropoietin
Purification: Antigen-affinity Purification
Homology: Sus scrofa (pig)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.