

Biotinylated Anti-rabbit IL-1 alpha Polyclonal Antibody

KPAB0202

Description

This Biotinylated Anti-rabbit IL-1 alpha Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0202
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Chicken
Reactivity: Rabbit
Immunogen: Recombinant Rabbit IL-1 alpha
Purification: Antigen-affinity Purification
Homology: Oryctolagus cuniculus (rabbit)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.