

Biotinylated Anti-swine IFN beta Polyclonal Antibody

KPAB0205

Description

This Biotinylated Anti-swine IFN beta Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0205
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Swine
Immunogen: Recombinant Porcine IFN beta
Purification: Antigen-affinity Purification
Homology: Phacochoerus africanus (Common warthog), Potamochoerus larvatus (bushpig), Potamochoerus porcus (red river hog), Sus scrofa (pig)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.